

<b>SonoAce 9900 Probe Guide</b>																
Transducer		C3-7IM	L4-7EV	L5-9ER	LS-12IM	P2-3AC	P2-5AC	P3-7AC	EC-4-9ES	EC-4-9IS	EC-4-9/10ED	C4-9/10ED	BCL4-9ES*	LIS-9EV	CL4-8EV	
Array Type		Curved	Linear	Linear	Linear	Phased	Phased	Phased	Curved	Curved	Curved	Curved	Curved/Linear	Linear	Curved	
Sector Angle/Field of View		77°	40mm	50mm	40mm	19mm	15mm	10mm	150°	150°	150°	150°	150°/55mm	40mm	65°	
Radius of Curvature		50mm	Flat	Flat	Flat	Flat	Flat	Flat	10mm	10mm	10mm	10mm	10mm/Flat	Flat	40mm	
Broadband Frequency Range		3MHz-7MHz	4MHz-7MHz	5MHz-9MHz	5MHz-12MHz	2MHz-3MHz	2MHz-5MHz	3MHz-7MHz	4MHz-9MHz	4MHz-9MHz	4MHz-9MHz	4MHz-9MHz	4MHz-9MHz	5MHz-9MHz	4MHz-9MHz	
<b>Application</b>																
Abdomen		●				●	●							●	●	
Gynecology		●							●	●	●					
OB		●							●	●	●					
Fetal Heart		●														
Urology							●		●	●	●		●			
Pediatric					●	●		●				●		●	●	
Small Part			●	●	●									●	●	
Breast			●	●	●											
Vascular			●	●	●		●									
Cardiology						●	●	●				●				
TCD						●										
Neonatal								●				●				
Intraoperative														●	●	
<b>Biopsy Guide</b>		Reusable	N/A	Reusable	Reusable	N/A	N/A	N/A	Reusable	Reusable	Reusable	N/A	Reusable	N/A	N/A	
		BP-C3-7M		VF-BP-PEC60	BP-LS-12M				BP-KIT-003	BP-KIT-003	BP-KIT-002-01		BP-BPL4-9EC			

<b>SonoAce 9900 Probe Guide</b>								
Transducer		VAW3-5	VAW4-7	VNW6-12	VDW5-8B	MPT4-7AL	2.0CW	4.0CW
Array Type		Curved	Curved	Linear	Curved	Phased	CW	CW
Sector Angle/Field of View		69°	69°	40mm	128°	10mm	N/A	N/A
Radius of Curvature		40mm	40mm	Flat	10mm	Flat	N/A	N/A
Broadband Frequency Range		3MHz-5MHz	4MHz-7MHz	6MHz-12MHz	5MHz-8MHz	4MHz-7MHz	2MHz	4MHz
<b>Application</b>								
Abdomen		●	●					
Gynecology		●	●		●			
OB		●	●		●			
Fetal Heart		●	●					
Urology					●			
Pediatric								●
Small Part				●				
Breast				●				
Vascular				●				
Cardiology						●	●	●
TCD								
Neonatal								
Intraoperative								
<b>Biopsy Guide</b>		Reusable	Reusable	Reusable	Reusable	N/A	N/A	N/A
				VF-BP-PEC54	VF-BP-PEC47			